



## STATION 2: Modified Bicron Analyst Measurements Form

Note: Only one exposed person should be in the room during measurements.

Date: \_\_\_\_\_ Time: \_\_\_\_\_ am / pm

Name of Exposed Individual: \_\_\_\_\_

1. Turn Modified Bicron Analyst on to x1 scale position.
2. Place the red dot on the survey meter probe against the exposed person's midsection.
3. Press the small round gray switch on the top right of the meter face. This button starts counting for one-minute intervals. A red light will come on during counting and will shut off at the end of each one-minute count.
4. Record the number displayed in the scaler box on the instrument meter face for each one-minute count.
5. Record the number shown in the scaler display box on Modified Bicron Analyst Measurements Form (R.P. Form # 97).
6. Collect and record at least 5 one-minute count measurements.

	Time of Day	Scaler Reading on Meter (counts per minute)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		

Measurements taken by: \_\_\_\_\_